

For Research Use Only

# TPRKB Monoclonal antibody, PBS Only

Catalog Number: 66402-1-PBS



## Basic Information

|   |  |   |
|---|--|---|
| <b>Catalog Number:</b><br>66402-1-PBS       | <b>GenBank Accession Number:</b><br>BC029492 | <b>Purification Method:</b><br>Protein A purification |
| <b>Size:</b><br>1mg/ml                      | <b>GeneID (NCBI):</b><br>51002               | <b>CloneNo.:</b><br>2D5B9                             |
| <b>Source:</b><br>Mouse                     | <b>UNIPROT ID:</b><br>Q9Y3C4                 |   |
| <b>Isotype:</b><br>IgG1                     | <b>Full Name:</b><br>TP53RK binding protein  |   |
| <b>Immunogen Catalog Number:</b><br>AG18770 | <b>Calculated MW:</b><br>175 aa, 20 kDa      |   |
|   | <b>Observed MW:</b><br>18-20 kDa             |   |

## Applications

**Tested Applications:**  
WB, IHC, IF/ICC, Indirect ELISA

**Species Specificity:**  
human

## Background Information

TPRKB, also named as CGI 121, belongs to the CGI121/TPRKB family. It interacts with TP53RK/PRPK and is widely expressed in all tissues.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

