

For Research Use Only

PDI Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66422-3-PBS



Basic Information

Catalog Number:

66422-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG1747

GenBank Accession Number:

BC014504

GeneID (NCBI):

5034

UNIPROT ID:

P07237

Full Name:

prolyl 4-hydroxylase, beta polypeptide

Calculated MW:

57 kDa

Purification Method:

Protein G purification

CloneNo.:

2E3B7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50798-1, PDI Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66422-2-PBS. Detection antibody: 66422-3-PBS. Standard:Ag1747. Range: 0.098-25 ng/mL



Cytometric bead array standard curve of MP50798-2, PDI Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66422-1-PBS. Detection antibody: 66422-3-PBS. Standard:Ag1747. Range: 0.098-1.563 ng/mL