

For Research Use Only

CD147 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66443-4-PBS



Basic Information

Catalog Number: 66443-4-PBS	GenBank Accession Number: BC009040	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 682	CloneNo.: 1E10H1
Source: Mouse	UNIPROT ID: P35613	
Isotype: IgG2b	Full Name: basigin (Ok blood group)	
Immunogen Catalog Number: AG2619	Calculated MW: 385 aa, 42 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50494-2, CD147 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66443-3-PBS. Detection antibody: 66443-4-PBS. Standard: Ag2619. Range: 0.195-25 ng/mL