

For Research Use Only

# IRF4 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 66451-2-PBS

## Basic Information

|   |  |   |
|---|--|---|
| <b>Catalog Number:</b><br>66451-2-PBS       | <b>GenBank Accession Number:</b><br>BC015752 | <b>Purification Method:</b><br>Protein G purification |
| <b>Size:</b><br>1mg/ml                      | <b>GeneID (NCBI):</b><br>3662                | <b>CloneNo.:</b><br>2B2C4                             |
| <b>Source:</b><br>Mouse                     | <b>ENSEMBL Gene ID:</b><br>ENSG00000137265   |   |
| <b>Isotype:</b><br>IgG1                     | <b>UNIPROT ID:</b><br>Q15306                 |   |
| <b>Immunogen Catalog Number:</b><br>AG25203 | <b>Full Name:</b><br>IRF 4                   |   |
|   | <b>Calculated MW:</b><br>52 kDa              |   |

## Applications

**Tested Applications:**  
Sandwich ELISA, Indirect ELISA

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.

**Storage Buffer:**  
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP50012-1, IRF4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66451-2-PBS. Detection antibody: HRP-conjugated 66451-1-PBS. Standard: Ag25203. Range: 0.02-20 ng/mL.



Sandwich ELISA standard curve of MP50012-2, IRF4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66451-2-PBS. Detection antibody: HRP-conjugated 66451-3-PBS. Standard: Ag25203. Range: 0.02-20 ng/mL.