

For Research Use Only

IRF4 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 66451-3-PBS

Basic Information

Catalog Number: 66451-3-PBS	GenBank Accession Number: BC015752	Purification Method: Protein G purification
Size: 1mg/ml	GeneID (NCBI): 3662	CloneNo.: 1F5C7
Source: Mouse	ENSEMBL Gene ID: ENSG00000137265	
Isotype: IgG1	UNIPROT ID: Q15306	
Immunogen Catalog Number: AG25203	Full Name: IRF 4	
	Calculated MW: 52 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50012-2, IRF4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66451-2-PBS. Detection antibody: HRP-conjugated 66451-3-PBS. Standard: Ag25203. Range: 0.02-20 ng/mL.