

For Research Use Only

# DCX Monoclonal antibody, PBS Only (Capture)



Catalog Number: 66462-2-PBS

## Basic Information

<b>Catalog Number:</b> 66462-2-PBS	<b>GenBank Accession Number:</b> BC027925	<b>Purification Method:</b> Protein G purification
<b>Size:</b> 1mg/ml	<b>GeneID (NCBI):</b> 1641	<b>CloneNo.:</b> 1F11B11
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> O43602	
<b>Isotype:</b> IgG1	<b>Full Name:</b> doublecortin	
<b>Immunogen Catalog Number:</b> AG4945	<b>Calculated MW:</b> 402 aa, 45 kDa	

## Applications

**Tested Applications:**  
Sandwich ELISA, Indirect ELISA

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP50026-1, DCX Monoclonal Matched Antibody Pair - BSA and Azide free. Capture antibody: 66462-2-PBS. Detection antibody: HRP-conjugated 66462-3-PBS. Standard: Ag4945. Range: 0.156-20 ng/mL