

For Research Use Only

# Rabbit Heavy chain of Rabbit IgG Monoclonal antibody, PBS Only



Catalog Number: 66467-1-PBS

## Basic Information

Catalog Number:

66467-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

GeneID (NCBI):

UNIPROT ID:

N/A

Full Name:

Observed MW:

50-55 kDa

Purification Method:

Protein G purification

CloneNo.:

1C9A5

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

rabbit

## Background Information

This antibody specifically recognizes the heavy chain of Rabbit IgG.

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



rabbit serum (10nL) were subjected to SDS PAGE followed by western blot with 66467-1-Ig (mouse McAb to Heavy chain of Rabbit IgG) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 66467-1-PBS in a different storage buffer formulation.