

For Research Use Only

# CASP3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66470-3-PBS



## Basic Information

Catalog Number: 66470-3-PBS	GenBank Accession Number: BC016926	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 836	CloneNo.: 1H9A7
Source: Mouse	UNIPROT ID: P42574	
Isotype: IgG2b	Full Name: caspase 3, apoptosis-related cysteine peptidase	
Immunogen Catalog Number: AG25029	Calculated MW: 277 aa, 32 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50501-1, CASP3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66470-3-PBS. Detection antibody: 66470-4-PBS. Standard: Ag25029. Range: 0.781-100 ng/mL.