

For Research Use Only

MEA1 Monoclonal antibody, PBS Only



Catalog Number: 66472-1-PBS

Basic Information

Catalog Number: 66472-1-PBS	GenBank Accession Number: BC001754	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 4201	CloneNo.: 1H8C2
Source: Mouse	UNIPROT ID: Q16626	
Isotype: IgG1	Full Name: male-enhanced antigen 1	
Immunogen Catalog Number: AG0375	Calculated MW: 20 kDa	
	Observed MW: 28 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, rat

Background Information

MEA1, also named as Male-enhanced antigen 1, is a 185 amino acid protein, which is highly expressed in testis. MEA1 may play an important role in spermatogenesis and/or testis development.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



rat testis tissue were subjected to SDS PAGE followed by western blot with 66472-1-Ig (MEA1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 66472-1-PBS in a different storage buffer formulation.