

For Research Use Only

# STAT2 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 66485-4-PBS

## Basic Information

|   |   |  |
|---|---|--|
| <b>Catalog Number:</b><br>66485-4-PBS       | <b>GenBank Accession Number:</b><br>BC051284                                    | <b>Purification Method:</b><br>Protein G Magarose purification |
| <b>Size:</b><br>1 mg/ml                     | <b>GeneID (NCBI):</b><br>6773   | <b>CloneNo.:</b><br>2D2B1                                      |
| <b>Source:</b><br>Mouse                     | <b>UNIPROT ID:</b><br>P52630  |  |
| <b>Isotype:</b><br>IgG1                     | <b>Full Name:</b><br>signal transducer and activator of transcription 2, 113kDa |  |
| <b>Immunogen Catalog Number:</b><br>AG10168 | <b>Calculated MW:</b><br>851 aa, 98 kDa   |  |

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50329-2, STAT2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66485-3-PBS. Detection antibody: 66485-4-PBS. Standard:Ag10168. Range: 1.563-200 ng/mL.