

For Research Use Only

Calretinin Monoclonal antibody, PBS Only (Detector)



Catalog Number: 66496-3-PBS

Basic Information

Catalog Number:

66496-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG2924

GenBank Accession Number:

BC015484

GeneID (NCBI):

794

UNIPROT ID:

P22676

Full Name:

calbindin 2

Calculated MW:

29 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2A1D5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50295-1, Calretinin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66496-2-PBS. Detection antibody: 66496-3-PBS. Standard: Ag2924. Range: 1.563-100 ng/mL.