

For Research Use Only

PCCB Monoclonal antibody, PBS Only (Capture)



Catalog Number: 66501-2-PBS

Basic Information

Catalog Number: 66501-2-PBS	GenBank Accession Number: BC013768	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 5096	CloneNo.: 1D6H12
Source: Mouse	UNIPROT ID: P05166	
Isotype: IgG2a	Full Name: propionyl Coenzyme A carboxylase, beta polypeptide	
Immunogen Catalog Number: AG1608	Calculated MW: 58 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50335-1, PCCB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66501-2-PBS. Detection antibody: 66501-3-PBS. Standard: Ag1608. Range: 1.563-100 ng/mL.



Cytometric bead array standard curve of MP50335-2, PCCB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66501-2-PBS. Detection antibody: 66501-1-PBS. Standard: Ag1608. Range: 0.781-100 ng/mL.