

For Research Use Only

PYCR1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66510-2-PBS



Basic Information

Catalog Number: 66510-2-PBS	GenBank Accession Number: BC022244	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 5831	CloneNo.: 1F4G1
Source: Mouse	UNIPROT ID: P32322	
Isotype: IgG1	Full Name: pyrroline-5-carboxylate reductase 1	
Immunogen Catalog Number: AG17919	Calculated MW: 319 aa, 33.8 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50786-1, PYCR1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66510-2-PBS. Detection antibody: 66510-1-PBS. Standard: Ag17919. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50786-3, PYCR1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66510-2-PBS. Detection antibody: 66510-4-PBS. Standard: Ag17919. Range: 0.098-100 ng/mL.