

For Research Use Only

PYCR1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66510-5-PBS



Basic Information

Catalog Number:

66510-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG17919

GenBank Accession Number:

BC022244

GeneID (NCBI):

5831

UNIPROT ID:

P32322

Full Name:

pyrroline-5-carboxylate reductase 1

Calculated MW:

319 aa, 33.8 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2A1C1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50786-4, PYCR1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66510-5-PBS. Detection antibody: 66510-4-PBS. Standard: Ag17919. Range: 0.098-100 ng/mL.