

For Research Use Only

OLIG2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66513-3-PBS



Basic Information

Catalog Number: 66513-3-PBS	GenBank Accession Number: BC047511	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 10215	CloneNo.: 1C2G3
Source: Mouse	UNIPROT ID: Q13516	
Isotype: IgG3	Full Name: oligodendrocyte lineage transcription factor 2	
Immunogen Catalog Number: AG18838	Calculated MW: 32 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51206-1, OLIG2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66513-2-PBS. Detection antibody: 66513-3-PBS. Standard: Ag18838. Range: 6.25-200 ng/mL