

For Research Use Only

Caspase 2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66517-2-PBS



Basic Information

Catalog Number:

66517-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG20141

GenBank Accession Number:

BC002427

GeneID (NCBI):

835

UNIPROT ID:

P42575

Full Name:

caspase 2, apoptosis-related cysteine peptidase

Calculated MW:

452 aa, 51 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

5B7F6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50435-1, Caspase 2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66517-2-PBS. Detection antibody: 66517-1-PBS. Standard: Ag20141. Range: 0.098-25 ng/mL