

For Research Use Only

SIAE Monoclonal antibody, PBS Only

Catalog Number: 66522-1-PBS



Basic Information

Catalog Number: 66522-1-PBS	GenBank Accession Number: BC068450	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 54414	CloneNo.: 1G5F7
Source: Mouse	UNIPROT ID: Q9HAT2	
Isotype: IgG1	Full Name: sialic acid acetyltransferase	
Immunogen Catalog Number: AG14624	Calculated MW: 58 kDa	
	Observed MW: 70 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human, mouse, rat, pig

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 66522-1-Ig (SIAE antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 66522-1-PBS in a different storage buffer formulation.



pig colon tissue were subjected to SDS PAGE followed by western blot with 66522-1-Ig (SIAE antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 66522-1-PBS in a different storage buffer formulation.