

For Research Use Only

IDO1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66528-3-PBS



Basic Information

Catalog Number: 66528-3-PBS	GenBank Accession Number: BC027882	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 3620	CloneNo.: 2B10B12
Source: Mouse	UNIPROT ID: P14902	
Isotype: IgG1	Full Name: indoleamine 2,3-dioxygenase 1	
Immunogen Catalog Number: AG3953	Calculated MW: 403 aa, 45 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50380-1, IDO1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66528-2-PBS. Detection antibody: 66528-3-PBS. Standard: Ag3953. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP50380-3, IDO1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66528-5-PBS. Detection antibody: 66528-3-PBS. Standard: Ag3953. Range: 6.25-100 ng/mL.