

For Research Use Only

MRPL45 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 66531-4-PBS



Basic Information

Catalog Number: 66531-4-PBS	GenBank Accession Number: BC006235	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 84311	CloneNo.: 1B9A11
Source: Mouse	UNIPROT ID: Q9BRJ2	
Isotype: IgG1	Full Name: mitochondrial ribosomal protein L45	
Immunogen Catalog Number: AG8207	Calculated MW: 35 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50680-2, MRPL45 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66531-4-PBS. Detection antibody: 66531-5-PBS. Standard: Ag8207. Range: 0.781-200 ng/mL



Cytometric bead array standard curve of MP50680-3, MRPL45 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66531-3-PBS. Detection antibody: 66531-4-PBS. Standard: Ag8207. Range: 3.125-200 ng/mL