

For Research Use Only

Cathepsin D Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66534-2-PBS



Basic Information

Catalog Number: 66534-2-PBS	GenBank Accession Number: BC016320	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 1509	CloneNo.: 2D6C10
Source: Mouse	UNIPROT ID: P07339	
Isotype: IgG2b	Full Name: cathepsin D	
Immunogen Catalog Number: AG15254	Calculated MW: 412 aa, 45 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51109-1, Cathepsin D Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66534-2-PBS. Detection antibody: 66534-3-PBS. Standard: Ag15254. Range: 0.781-100 ng/mL.