

For Research Use Only

Calponin Monoclonal antibody, PBS Only (Capture)



Catalog Number: 66540-2-PBS

Basic Information

Catalog Number:

66540-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG21384

GenBank Accession Number:

BC022015

GeneID (NCBI):

1264

UNIPROT ID:

P51911

Full Name:

calponin 1, basic, smooth muscle

Calculated MW:

33 kDa

Purification Method:

Protein G purification

CloneNo.:

1B3C9

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50017-1, Calponin Monoclonal Matched Antibody Pair - BSA and Azide free. Capture antibody: 66540-2-PBS. Detection antibody: HRP-conjugated 66540-1-PBS. Standard: Ag21384. Range: 0.313-20 ng/mL.