

For Research Use Only

TOP2A Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66541-3-PBS



Basic Information

Catalog Number:

66541-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG19551

GenBank Accession Number:

BC140791

GeneID (NCBI):

7153

UNIPROT ID:

P11388

Full Name:

topoisomerase (DNA) II alpha 170kDa

Calculated MW:

1531 aa, 174 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

4A3A12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51073-1, TOP2A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66541-2-PBS. Detection antibody: 66541-3-PBS. Standard:Ag19551. Range: 0.195-100 ng/mL



Cytometric bead array standard curve of MP51073-2, TOP2A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66541-4-PBS. Detection antibody: 66541-3-PBS. Standard:Ag19551. Range: 0.195-100 ng/mL