

For Research Use Only

# KI67 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 66555-5-PBS

## Basic Information

Catalog Number: 66555-5-PBS	GenBank Accession Number: NM_002417	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 4288	CloneNo.: 1D8E4
Source: Mouse	UNIPROT ID: P46013	
Isotype: IgG1	Full Name: antigen identified by monoclonal antibody Ki-67	
Immunogen Catalog Number: AG26266	Calculated MW: 359 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50237-3, KI67 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66555-5-PBS. Detection antibody: 66555-6-PBS. Standard: Ag26266. Range: 0.098-50 ng/mL