

For Research Use Only

HE4 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66557-2-PBS



Basic Information

Catalog Number:

66557-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG5936

GenBank Accession Number:

BC046106

GeneID (NCBI):

10406

UNIPROT ID:

Q14508

Full Name:

WAP four-disulfide core domain 2

Calculated MW:

13 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

5D4E5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51083-1, HE4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66557-2-PBS. Detection antibody: 66557-3-PBS. Standard:Ag5936. Range: 0.195-100 ng/mL.