

For Research Use Only

HE4 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66557-4-PBS



Basic Information

Catalog Number: 66557-4-PBS	GenBank Accession Number: BC046106	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 10406	CloneNo.: 4B9D3
Source: Mouse	UNIPROT ID: Q14508	
Isotype: IgG2a	Full Name: WAP four-disulfide core domain 2	
Immunogen Catalog Number: AG5936	Calculated MW: 13 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51083-3, HE4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66557-4-PBS. Detection antibody: 66557-5-PBS. Standard: Ag5936. Range: 0.781-100 ng/mL.