

For Research Use Only

SMARCA4 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66561-2-PBS



Basic Information

Catalog Number:

66561-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG16256

GenBank Accession Number:

BC150298

GeneID (NCBI):

6597

UNIPROT ID:

P51532

Full Name:

SWI/SNF related, matrix associated, actin dependent regulator of chromatin, subfamily a, member 4

Calculated MW:

1647 aa, 185 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2E7B8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51084-1, SMARCA4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66561-2-PBS. Detection antibody: 66561-3-PBS. Standard: Ag16256. Range: 0.098-100 ng/mL.