

For Research Use Only

ATF6 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66563-3-PBS



Basic Information

Catalog Number:

66563-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG21456

GenBank Accession Number:

BC014969

GeneID (NCBI):

22926

UNIPROT ID:

P18850

Full Name:

activating transcription factor 6

Calculated MW:

75 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

3A9B12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51101-1, ATF6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66563-2-PBS. Detection antibody: 66563-3-PBS. Standard:ATF6. Range: 0.098-50 ng/mL