

For Research Use Only

CD151 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 66567-3-PBS

Basic Information

Catalog Number: 66567-3-PBS	GenBank Accession Number: BC001374	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 977	CloneNo.: 2F10E1
Source: Mouse	UNIPROT ID: P48509	
Isotype: IgG1	Full Name: CD151 molecule (Raph blood group)	
Immunogen Catalog Number: AG25985	Calculated MW: 28 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50304-2, CD151 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66567-3-PBS. Detection antibody: 66567-4-PBS. Standard: Ag25985. Range: 6.25-200 ng/mL