

For Research Use Only

GRP78,BIP Monoclonal antibody, PBS Only (Capture)

Catalog Number:66574-2-PBS



Basic Information

Catalog Number:

66574-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG2188

GenBank Accession Number:

BC020235

GeneID (NCBI):

3309

UNIPROT ID:

P11021

Full Name:

heat shock 70kDa protein 5 (glucose-regulated protein, 78kDa)

Calculated MW:

78 kDa

Purification Method:

Protein G purification

CloneNo.:

1G1F4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50855-1, GRP78,BIP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66574-2-PBS. Detection antibody: 66574-3-PBS. Standard:Ag2188. Range: 1.563-25 ng/mL.