

For Research Use Only

GRP78,BIP Monoclonal antibody, PBS Only (Detector)

Catalog Number:66574-3-PBS



Basic Information

Catalog Number: 66574-3-PBS	GenBank Accession Number: BC020235	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 3309	CloneNo.: 1E5E6
Source: Mouse	UNIPROT ID: P11021	
Isotype: IgG1	Full Name: heat shock 70kDa protein 5 (glucose-regulated protein, 78kDa)	
Immunogen Catalog Number: AG2188	Calculated MW: 78 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50855-1, GRP78,BIP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66574-2-PBS. Detection antibody: 66574-3-PBS. Standard:Ag2188. Range: 1.563-25 ng/mL.