

For Research Use Only

TOM70 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66593-4-PBS



Basic Information

Catalog Number: 66593-4-PBS	GenBank Accession Number: BC052994	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 9868	CloneNo.: 2A2H10
Source: Mouse	UNIPROT ID: O94826	
Isotype: IgG2a	Full Name: translocase of outer mitochondrial membrane 70 homolog A (S. cerevisiae)	
Immunogen Catalog Number: AG6027	Calculated MW: 67 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50612-2, TOM70 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66593-2-PBS. Detection antibody: 66593-4-PBS. Standard: Ag6027. Range: 1.563-100 ng/mL.