

For Research Use Only

# MRPL23 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66595-2-PBS



## Basic Information

Catalog Number: 66595-2-PBS	GenBank Accession Number: BC027710	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 6150	CloneNo.: 3C3C3
Source: Mouse	UNIPROT ID: Q16540	
Isotype: IgG2a	Full Name: mitochondrial ribosomal protein L23	
Immunogen Catalog Number: AG2220	Calculated MW: 153 aa, 18 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50460-1, MRPL23 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66595-2-PBS. Detection antibody: 66595-3-PBS. Standard: Ag2220. Range: 0.391-100 ng/mL.