

For Research Use Only

# MRPL23 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66595-3-PBS



## Basic Information

|                                     |   |  |
|-------------------------------------|---|--|
| Catalog Number:<br>66595-3-PBS      | GenBank Accession Number:<br>BC027710             | Purification Method:<br>Protein G purification |
| Size:<br>1 mg/ml                    | GeneID (NCBI):<br>6150                            | CloneNo.:<br>3E2C4                             |
| Source:<br>Mouse                    | UNIPROT ID:<br>Q16540                             |  |
| Isotype:<br>IgG1                    | Full Name:<br>mitochondrial ribosomal protein L23 |  |
| Immunogen Catalog Number:<br>AG2220 | Calculated MW:<br>153 aa, 18 kDa                  |  |

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50460-1, MRPL23 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66595-2-PBS. Detection antibody: 66595-3-PBS. Standard:Ag2220. Range: 0.391-100 ng/mL.



Cytometric bead array standard curve of MP50460-2, MRPL23 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66595-4-PBS. Detection antibody: 66595-3-PBS. Standard:Ag2220. Range: 0.391-100 ng/mL.