

For Research Use Only

# CHCHD6 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66597-2-PBS



## Basic Information

Catalog Number:

66597-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG14700

GenBank Accession Number:

BC006123

GeneID (NCBI):

84303

UNIPROT ID:

Q9BRQ6

Full Name:

coiled-coil-helix-coiled-coil-helix domain containing 6

Calculated MW:

235 aa, 26 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2A10B5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50477-1, CHCHD6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66597-2-PBS. Detection antibody: 66597-3-PBS. Standard: Ag14700. Range: 0.391-100 ng/mL.



Cytometric bead array standard curve of MP50477-3, CHCHD6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66597-2-PBS. Detection antibody: 66597-4-PBS. Standard: Ag14700. Range: 0.781-100 ng/mL.