

For Research Use Only

CHCHD6 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66597-4-PBS



Basic Information

Catalog Number:

66597-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG14700

GenBank Accession Number:

BC006123

GeneID (NCBI):

84303

UNIPROT ID:

Q9BRQ6

Full Name:

coiled-coil-helix-coiled-coil-helix domain containing 6

Calculated MW:

235 aa, 26 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1B6D11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50477-2, CHCHD6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66597-1-PBS. Detection antibody: 66597-4-PBS. Standard: Ag14700. Range: 1.563-100 ng/mL.



Cytometric bead array standard curve of MP50477-3, CHCHD6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66597-2-PBS. Detection antibody: 66597-4-PBS. Standard: Ag14700. Range: 0.781-100 ng/mL.