

For Research Use Only

# ATPB Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66600-2-PBS



## Basic Information

<b>Catalog Number:</b> 66600-2-PBS	<b>GenBank Accession Number:</b> BC016512	<b>Purification Method:</b> Protein G Magarose purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 506	<b>CloneNo.:</b> 1G3A4
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> P06576	
<b>Isotype:</b> IgG1	<b>Full Name:</b> ATP synthase, H <sup>+</sup> transporting, mitochondrial F1 complex, beta polypeptide	
<b>Immunogen Catalog Number:</b> AG11177	<b>Calculated MW:</b> 529 aa, 57 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage


**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:  
T: 4006900926      E: Proteintech-CN@ptglab.com      W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50492-1, ATPB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66600-2-PBS. Detection antibody: 66600-3-PBS. Standard: Ag11177. Range: 0.195-3.125 ng/mL.