

For Research Use Only

ATPB Monoclonal antibody, PBS Only (Capture)

Catalog Number:66600-4-PBS



Basic Information

Catalog Number: 66600-4-PBS	GenBank Accession Number: BC016512	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 506	CloneNo.: 2F5A11
Source: Mouse	UNIPROT ID: P06576	
Isotype: IgG1	Full Name: ATP synthase, H ⁺ transporting, mitochondrial F1 complex, beta polypeptide	
Immunogen Catalog Number: AG11177	Calculated MW: 529 aa, 57 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50492-2, ATPB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66600-4-PBS. Detection antibody: 66600-5-PBS. Standard: Ag11177. Range: 0.195-3.125 ng/mL.