

For Research Use Only

ATPB Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66600-7-PBS



Basic Information

Catalog Number:

66600-7-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG11177

GenBank Accession Number:

BC016512

GeneID (NCBI):

506

UNIPROT ID:

P06576

Full Name:

ATP synthase, H⁺ transporting,
mitochondrial F1 complex, beta
polypeptide

Calculated MW:

529 aa, 57 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3E2A9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50492-3, ATPB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66600-6-PBS. Detection antibody: 66600-7-PBS. Standard: Ag11177. Range: 1.563-100 ng/mL.