

For Research Use Only

ERMN Monoclonal antibody, PBS Only



Catalog Number: 66605-1-PBS

Basic Information

Catalog Number: 66605-1-PBS	GenBank Accession Number: BC026345	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 57471	CloneNo.: 2B1C1
Source: Mouse	UNIPROT ID: Q8TAM6	
Isotype: IgG2b	Full Name: ermin, ERM-like protein	
Immunogen Catalog Number: AG27040	Observed MW: 50 kDa	

Applications

Tested Applications:
WB, Indirect ELISA, IHC

Species Specificity:
Human, mouse, rat, pig

Background Information

ERMN is associated with cytoskeletal rearrangements and the stability of myelin sheath. A study has shown that ERMN is expressed specifically in oligodendrocytes in the central nervous system (CNS), and is involved in the formation of myelin (PMID: 20934411). The predicted MW of ERMN is 33 kDa, while western blot analyses of different tissue lysates using 66605-1-Ig antibody to ERMN detected a 50 kDa single band protein as a result of phosphorylation modification.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various tissues were subjected to SDS PAGE followed by western blot with 66605-1-Ig (ERMN antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 66605-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 66605-1-Ig (ERMN antibody) at dilution of 1:300 (under 10x lens. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 66605-1-PBS in a different storage buffer formulation.