

For Research Use Only

TNFSF11/RANKL Monoclonal antibody, PBS Only (Detector)

Catalog Number:66610-2-PBS



Basic Information

Catalog Number: 66610-2-PBS	GenBank Accession Number: BC074890	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 8600	CloneNo.: 1H4A2
Source: Mouse	UNIPROT ID: O14788	
Isotype: IgG1	Full Name: tumor necrosis factor (ligand) superfamily, member 11	
Immunogen Catalog Number: AG19975	Calculated MW: 317 aa, 35 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50526-1, RANKL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66610-1-PBS. Detection antibody: 66610-2-PBS. Standard: Ag19975. Range: 1.563-100 ng/mL.