

For Research Use Only

GRIK2 Monoclonal antibody, PBS Only



Catalog Number: 66631-2-PBS

Basic Information

Catalog Number: 66631-2-PBS	GenBank Accession Number: BC037954	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 2898	CloneNo.: 1E5G6
Source: Mouse	UNIPROT ID: Q13002	
Isotype: IgG1	Full Name: glutamate receptor, ionotropic, kainate 2	
Immunogen Catalog Number: AG25308	Calculated MW: 908 aa, 103 kDa	
	Observed MW: 103-115 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human, Mouse, Rat, Pig

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 66631-2-Ig (GRIK2 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 66631-2-PBS in a different storage buffer formulation.