

For Research Use Only

# p70(S6K) Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66638-2-PBS



## Basic Information

Catalog Number:

66638-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG5883

GenBank Accession Number:

BC053365

GeneID (NCBI):

6198

UNIPROT ID:

P23443

Full Name:

ribosomal protein S6 kinase, 70kDa, polypeptide 1

Calculated MW:

59 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1E9A5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50673-1, p70(S6K) Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66638-1-PBS. Detection antibody: 66638-2-PBS. Standard: Ag5883. Range: 1.563-100 ng/mL.