

For Research Use Only

CD302 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66640-2-PBS



Basic Information

Catalog Number: 66640-2-PBS	GenBank Accession Number: BC020646	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 9936	CloneNo.: 1C2H12
Source: Mouse	UNIPROT ID: Q8IX05	
Isotype: IgG1	Full Name: CD302 molecule	
Immunogen Catalog Number: AG9873	Calculated MW: 170aa, 19 kDa; 232aa, 26 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50401-1, CD302 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66640-2-PBS. Detection antibody: 66640-3-PBS. Standard: Ag9873. Range: 3.125-100 ng/mL.