

For Research Use Only

TLR2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66645-2-PBS



Basic Information

Catalog Number: 66645-2-PBS	GenBank Accession Number: BC033756	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 7097	CloneNo.: 2F4D5
Source: Mouse	UNIPROT ID: O60603	
Isotype: IgG1	Full Name: toll-like receptor 2	
Immunogen Catalog Number: AG25626	Calculated MW: 784 aa, 90 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50388-1, TLR2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66645-1-PBS. Detection antibody: 66645-2-PBS. Standard: Ag25626. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP50388-2, TLR2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66645-3-PBS. Detection antibody: 66645-2-PBS. Standard: Ag25626. Range: 0.781-100 ng/mL.