

For Research Use Only

# POU2AF1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66659-3-PBS



## Basic Information

Catalog Number:

66659-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG23613

GenBank Accession Number:

BC032549

GeneID (NCBI):

5450

ENSEMBL Gene ID:

ENSG00000110777

UNIPROT ID:

Q16633

Full Name:

POU class 2 associating factor 1

Calculated MW:

256 aa, 27 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

4G3A2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP51058-2, POU2AF1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66659-1-PBS. Detection antibody: 66659-3-PBS. Standard:Ag23613. Range: 0.781-200 ng/mL