

For Research Use Only

DDX54 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66664-2-PBS



Basic Information

Catalog Number:

66664-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG25289

GenBank Accession Number:

BC156669

GeneID (NCBI):

79039

UNIPROT ID:

Q8TDD1

Full Name:

DEAD (Asp-Glu-Ala-Asp) box
polypeptide 54

Calculated MW:

98 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2F6B2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51115-1, DDX54 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66664-2-PBS. Detection antibody: 66664-3-PBS. Standard: Ag25289. Range: 3.125-50 ng/mL