

For Research Use Only

# PDIA6 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66669-3-PBS



## Basic Information

Catalog Number:

66669-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG27608

GenBank Accession Number:

BC001312

GeneID (NCBI):

10130

UNIPROT ID:

Q15084

Full Name:

protein disulfide isomerase family A, member 6

Calculated MW:

48 kDa, 54 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1A9C7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50818-1, PDIA6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66669-2-PBS. Detection antibody: 66669-3-PBS. Standard: Ag27608. Range: 0.781-100 ng/mL.