

For Research Use Only

PDIA6 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66669-5-PBS



Basic Information

Catalog Number:

66669-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG27608

GenBank Accession Number:

BC001312

GeneID (NCBI):

10130

UNIPROT ID:

Q15084

Full Name:

protein disulfide isomerase family A, member 6

Calculated MW:

48 kDa, 54 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2C2C8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50818-2, PDIA6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66669-4-PBS. Detection antibody: 66669-5-PBS. Standard: Ag27608. Range: 3.125-100 ng/mL.



Cytometric bead array standard curve of MP50818-3, PDIA6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66669-6-PBS. Detection antibody: 66669-5-PBS. Standard: Ag27608. Range: 6.25-100 ng/mL.