

For Research Use Only

MLKL Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66675-4-PBS



Basic Information

Catalog Number: 66675-4-PBS	GenBank Accession Number: BC028141	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 197259	CloneNo.: 3C2A5
Source: Mouse	UNIPROT ID: Q8NB16	
Isotype: IgG1	Full Name: mixed lineage kinase domain-like	
Immunogen Catalog Number: AG15188	Calculated MW: 471 aa, 54 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51035-2, MLKL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66675-2-PBS. Detection antibody: 66675-4-PBS. Standard: Ag15188. Range: 3.125-100 ng/mL.