

For Research Use Only

AKR7A2 Monoclonal antibody, PBS Only (Detector)

Catalog Number:66677-4-PBS



Basic Information

Catalog Number:

66677-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG26105

GenBank Accession Number:

BC007352

GeneID (NCBI):

8574

UNIPROT ID:

O43488

Full Name:

aldo-keto reductase family 7,
member A2 (aflatoxin aldehyde
reductase)

Calculated MW:

40 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2F1A9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50657-2, AKR7A2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66677-2-PBS. Detection antibody: 66677-4-PBS. Standard: Ag26105. Range: 1.563-100 ng/mL.